

Kinase Assay Using Radioisotope (Autoradiography)

Preparation of assay conditions, a 20µL reaction containing the following:

20 mM Hepes KOH (pH7.5)

10 mM MgCl₂

1 mM DTT

100 µM [gamma 32P] ATP (0.5µCi)

4 µL of Abnova's enzyme and substrate cocktail (2 µg each of histone H1, histone H3, MBP and beta casein)

Step-by-step procedure:

1. Start the reaction by adding 4 µL of [gamma 32P] ATP (0.5µCi)
2. Incubate for 60 minutes at 30°C.
3. Terminate the reaction by adding 7 µL of 4X sample buffer.
4. Subsequently boil the sample.
5. Analyze the 15 µL of sample on SDS-PAGE gel.
6. Stain, destain, and then dry the gel.
7. Take an autoradiograph by exposing the gel for 24 hours.

Kinase substrates:

H3 histone, family 3A - **calf thymus**

Molecular weight: approx. 16.5-17 kDa on SDS-PAGE

Myelin basic protein - **bovine brain**

Molecular weight: approx. 21-22 kDa on SDS-PAGE

Histone H1.0 - **calf thymus**

Molecular weight: approx. 32-33 kDa on SDS-PAGE

Beta-casein protein- **bovine milk**

Molecular weight: approx. 33-34 kDa on SDS-PAGE

Kinase Assay Using Radioisotope (Autoradiography)

Preparation of assay conditions, a 20µL reaction containing the following:

20 mM Hepes KOH (pH7.5)

10 mM MgCl₂

1 mM DTT

0.1 mM EGTA

2.5 mM CaCl₂

250 ng/ml calmodulin

100 µM [gamma 32P] ATP (0.5µCi)

4 µL of Abnova's enzyme and substrate cocktail (2 µg each of histone H1, histone H3, MBP and beta casein)

Step-by-step procedure:

1. Start the reaction by adding 4 µL of [gamma 32P] ATP (0.5µCi)
2. Incubate for 60 minutes at 30°C.



Abnova Products are distributed in the UK and Ireland by Caltag Medsystems Ltd,
Botolph Claydon, Buckingham, MK18 2LR, UK.

Tel: Freephone (UK) 0800 279 9113, or +44 (0) 1296 715459

Email: office@caltagmedsystems.co.uk Web: www.caltagmedsystems.co.uk

3. Terminate the reaction by adding 7 μ L of 4X sample buffer.
4. Subsequently boil the sample.
5. Analyze the 15 μ L of sample on SDS-PAGE gel.
6. Stain, destain, and then dry the gel.
7. Take an autoradiograph by exposing the gel for 24 hours.

Kinase substrates:

Beta-casein protein- **bovine milk**

Molecular weight: approx. 33-34 kDa on SDS-PAGE

Myelin basic protein - **bovine brain**

Molecular weight: approx. 21-22 kDa on SDS-PAGE

Histone H1.0 - **calf thymus**

Molecular weight: approx. 32-33 kDa on SDS-PAGE

H3 histone, family 3A - **calf thymus**

Molecular weight: approx. 16.5-17 kDa on SDS-PAGE

Phosphatase Assay Using Fluorescent Substrates

Preparation of assay conditions, a 50 μ L reaction containing the following:

100 mM Tris HCl (pH 7.5)*

40 mM NaCl

1 mM DTT

0.1 mM OMFP (3-o-methylfluorescein phosphate)

5 μ L of Abnova's enzyme in each well of a 96-well microplate

Step-by-step procedure:

1. Measure fluorescence intensity of sample for 60 minutes at 5-minute intervals using a microtiter plate fluorometer with excitation at 485 nm and emission at 535 nm at 30°C.

2. Measure and calculate the rate of reaction while the velocity of the reaction remains constant.

*50 mM Mops KOH pH 6.0 for acid phosphatase.

Equipment(s):

Microtiter plate fluorometer (Arvo HTS; Perkin Elmer)



Abnova Products are distributed in the UK and Ireland by Caltag Medsystems Ltd,
Botolph Claydon, Buckingham, MK18 2LR, UK.

Tel: Freephone (UK) 0800 279 9113, or +44 (0) 1296 715459

Email: office@caltagmedsystems.co.uk Web: www.caltagmedsystems.co.uk