

## Recommended Immunohistochemistry Protocol:

The following protocol is a recommendation only, and MBL makes no guarantee of the results:

- Tissue Preparation: Formalin fixation and embedding in paraffin wax
- Tissue Sectioning: Make 4- $\mu$ m sections and place on pre-cleaned and charged microscope slides. Heat in a tissue-drying oven for 45 minutes at 60°C.
- Deparaffinization: Wash dry slides in 3 changes of xylene – 5 minutes each @ RT
- Rehydration: Wash slides in 3 changes of 100% alcohol – 3 minutes each @ RT  
Wash slides in 2 changes of 95% alcohol – 3 minutes each @ RT  
Wash slides in 1 change of 80% alcohol – 3 minutes @ RT  
Rinse slides in gentle running distilled water – 5 minutes @ RT
- Antigen retrieval: Steam slides in 0.01 M sodium citrate buffer, pH 6.0 at 99-100°C - 20 minutes  
Remove from heat and let stand at room temperature in buffer - 20 minutes  
Rinse in 1X TBS with Tween (TBST) – 1 minute @ RT
- Immunostaining: (Do not allow tissues to dry at any time during the staining procedure)  
Apply a universal protein block – 20 minutes @ RT  
Drain protein block from slides, apply diluted primary antibody – 45 minutes @ RT  
Rinse slides in 1X TBST - 1 minute @ RT  
Apply a biotinylated anti-rabbit IgG (H+L) secondary – 30 minutes @ RT  
Rinse slides in 1X TBST - 1 minute @ RT  
Apply alkaline phosphatase streptavidin – 30 minutes @ RT  
Rinse slides in 1X TBST - 1 minute @ RT  
Apply alkaline phosphatase chromogen substrate – 30 minutes @ RT  
Wash slides in distilled water – 1 minute @ RT
- Dehydrate: (This method should only be used if the chromogen substrate is alcohol insoluble (e.g. Vector Red, DAB))  
Wash slides in 2 changes of 80% alcohol – 1 minute each @ RT  
Wash slides in 2 changes of 95% alcohol – 1 minute each @ RT  
Wash slides in 3 changes of 100% alcohol – 1 minute each @ RT  
Wash slides in 3 changes of xylene – 1 minute each @ RT  
Apply coverslip

- The sodium citrate buffer, TBST, universal block and antibody diluent are purchased from Dako, and the secondary antibodies and chromogen (usually Vector Red) are obtained from Vector.
- Our tissues are fixed in 10% neutral buffered formalin. The ratio of buffered formalin to tissue volume is 10:1 for approximately 18 hours. Formalin has a very slow diffusion coefficient so the tissue needs to be no more than about the thickness of a nickel.